

Mechanical Lift Checklist



Lift Number _____ Month _____ Year _____ (Check One) Sit to Stand -OR- Full Body/Hoyer

Area Check	Week 1	Week 2	Week 3	Week 4	Week 5
Hoist Check Raise and Lower Function.					
Wheels Move Freely. In good overall condition.					
Cradle Assembly and Base In good condition. Adequate width movement.					
Carrier and Boom Attachment Point In good condition.					
Front Legs In good condition.					
Nuts, Bolts, Weld Points In good condition. Connections are tight.					
Overall No cracks, damage, or sharp surfaces.					
Functionality Full function test completed. Functioning well.					
Initials					
Date					